

MODIFICATIONS

Problem	Modification(s)
Skin irritation or rubbing	<ol style="list-style-type: none"> 1. Wear toe socks. 2. Ensure footwear is wide enough 3. Use one of the following modifications to trim away bothersome material
Too much spread on pinky toe	<ol style="list-style-type: none"> 1. Cut off half of the pinky spacer. 2. Cut off the full pinky spacers. Shown in GREEN
Correct Toes slide off foot	<ol style="list-style-type: none"> 1. Try toe socks, or conventional sock over CT (stretch socks' toe seams first). 2. Trim big toe and/or pinky toe pillar, to allow CT to slide further on. Shown in BLUE
Too tight on toe(s)	<ol style="list-style-type: none"> 1. Try a larger size of CT. 2. Cut a small hole on the top and/or bottom. Shown in PURPLE
Bunion Tailor's Bunion	<p>Add a shim in the appropriate CT slot. See YELLOW arrow.</p> <p><i>*Only add shims after 3-4 months of regular Correct Toes use.</i></p> <ul style="list-style-type: none"> • <i>What's a shim? Take an old shoe liner and cut a piece to fit the slot in your Correct Toes.</i>
Hammertoe	<ol style="list-style-type: none"> 1. Cut a small hole in CT over your hammertoes. Shown in PURPLE 2. Trim away the interspace. Shown in MAGENTA <p><i>*Correct Toes are only effective for hammertoes that are flexible (can be manually straightened).</i></p>
Corn	<ol style="list-style-type: none"> 1. Apply Moleskin 2. Cut a small hole in the top of CT over your affected toe. Shown in PURPLE or MAGENTA 3. Trim away a portion of the pillar. Shown in red
Webbedtoes	Trim away a portion, or entirely remove the pillar where it conflicts with webbed toes. Shown in ORANGE

